




**FRESENIUS  
KABI**

caring for life



## Calculate and Prepare From Anywhere With Our New RBCX Estimator App

**Experience our new estimator  
app designed to support Amicus®  
Apheresis Operators**

- Provides a mobile solution for estimating the volume of replacement fluid needed for each patient in advance of the RBCX procedure
- Easily estimate the FCR that can be achieved using Amicus to perform the RBCX procedure<sup>1</sup>
- Access the app on the go without the need for a network connection or the need to be in front of the Amicus device

<sup>1</sup>FCR estimation is based on the physician prescribed replacement volumes

### Download the RBCX Estimator App Today!

Visit [gentlescience.fresenius-kabi.us](http://gentlescience.fresenius-kabi.us) and download the app today. Use your iPhone, iPad, or Android device.



### Resources at Your Fingertips - Anytime, anywhere, & From Any Mobile Device.

Our new app provides you with access to estimated volumes anywhere, anytime, without the need to be in front of the device.

### Questions about our RBCX Estimator App?

Ask your Fresenius Kabi Clinical Specialist for more information, or call our Clinical Support Line at 800-333-6925, or email us at [GentleScience@Fresenius-Kabi.com](mailto:GentleScience@Fresenius-Kabi.com)

Please see reverse side for limitations regarding the use of this App

This tool is not intended to be used to make medical decisions. It is intended as an aid for obtaining an approximation of the Amicus Separator calculations. The patient and prescription parameters must still be entered on the Amicus Separator at the time of the procedure for actual calculation of replacement volumes.

This app is for exclusive use in estimating volumes for an Amicus Separator RBCX procedure. This app is not to be used to estimate volumes for other devices used for RBCX procedures.

Refer to the Amicus Separator Operator's Manual for the full list of warnings and precautions associated with the use of the Amicus device.

Amicus is a registered trademark of Fresenius Kabi AG and/or its affiliates.



**FRESENIUS  
KABI**  
caring for life

Fresenius Kabi USA, LLC  
Three Corporate Drive  
Lake Zurich, Illinois 60047  
800-333-6925

Copyright © 2019 Fresenius Kabi AG  
0001-AMTH-01-06/19